


HEARING PRODUCTS FOR A BETTER LIFESTYLE



TOP PRODUCTS









ABOUT SARABEC

SARABEC® HAS OVER 40 YEARS **OF EXPERIENCE WORKING WITH PEOPLE WHO HAVE A HEARING IMPAIRMENT.**

At Sarabec, we aim to make life a little easier for those with hearing difficulties or hearing loss by offering simple hearing solutions that can be used with or without hearing aids. We provide high-quality, reliable assistive devices, products, and accessories to help individuals stay connected with friends and relatives, know what's happening at home, and enjoy recreational time.

Our products include listening assistance for TV and radio. So, while those around you listen to a volume of their choice, you can select a volume via the listening device that only you can hear. Our telephones come in a range of amplified corded, cordless, and mobile options, and our 'alerting' systems like extra loud doorbells, alarms, alarm clocks, and baby monitors aim to help keep you safe and independent.

At Sarabec, we offer a simple and easy buying process when you are ordering online or over the phone. With unlimited technical support, we are always on hand to answer any questions or queries you may have, whether it's before or after making a purchase.

1 MONTH TRIAL

With our 1-month trial policy, you can be sure that the products you choose are right for you by trying them in the comfort of your own home. If you are not happy with the product, simply contact us and we will let you know the best way to return the product to obtain a refund. We cannot give you a

refund if goods, when returned, are found to be damaged, incomplete or have been tampered with.

VAT RELIEF

If you have a hearing loss and are buying for your own use, you can obtain VAT Relief on certain products in this catalogue. Products that are eligible for VAT Relief are indicated as shown below.

Note: Some accessories if not bought as part of a system cannot be given VAT Relief.



This product qualifies for VAT Relief.

Relief.


This product does not qualify for VAT

To obtain VAT Relief you must complete a VAT Relief Certificate at the time of your order. When ordering online please tick the relevant section. VAT refunds cannot be given afterwards. If you are unsure, or want to know more, please contact us.

EVERYDAY LIFE WITH A PACEMAKER

Pacemaker manufacturers give detailed instructions about possible sources of interference from products that you should avoid. Please refer to your booklet about your type of pacemaker. If you have any further questions about what you can and can't do, please contact your pacemaker clinic. Further information is given on our website.

Symbols used throughout this catalogue:





CONTENTS

01	TV & AUDIO
03	PERSONAL LISTENERS
04	ALERTING - SIGNOLUX
06	ALERTING - BABY MONITOR
07	ALERTING - CHIMEFLASH
09	ALERTING - ALARM CLOCKS
10	PHONES
13	TINNITUS SNOOZEBANDS
14	HOW TO CONNECT TO YOUR TV
16	ORDERING INFORMATION

TV & Audio

Sarabec's TV and audio systems aim to make life easier for anyone who has difficulty hearing their television, radio, Hi-Fi or any other audio.

Amplify the sounds you want to hear whilst reducing disturbing background sounds. For use with or without a hearing aid, these systems bring sounds direct to your ears without disturbing others.

For information on connecting your devices to your TV, head to page 14.

Sonumaxx 2.4 PR Headset system



Transmitter, pocket receiver and plug-in stethoscope headset

Sonumaxx 2.4	PR Headset system
Order Code: 230-440	Colour: Black/Grey
Price: (Ex. VAT) £119.90 Price: (Inc. VAT) £143.88	Size : Base: 125mm Ø x 40mm Receiver : 225 x 125 x 25mm
VAT Relief Ves	

Sonumaxx 2.4 PR Neckloop System

Transmitter, pocket receiver and plug-in neckloop

Sonumaxx 2.4 PR Neckloop System

Order Code: 230-441	Colour: Black/Grey
Price: (Ex. VAT) £119.90	Size : Base: 125mm Ø x 40mm
Price: (Inc. VAT) £143.88	Receiver : 225 x 125 x 25mm







The LA240 Room Loop and Loop Pad System provide a practical solution for hearing aid users to more easily listen to their TV or Audio equipment via the "T" or "Loop" program of their hearing aid.

With no direct connection between the user and a TV or Audio equipment, the user is able to listen comfortably to a TV or audio equipment without the need to use a separate receiver.

The **LA240 Loop System** picks up the sound from a TV via a direct connection or microphone and feeds it to the amplifier. The sound is then passed to the **LA240 Loop Pad** which in turn transmits the sound inductively, "magnetic sound waves", to a hearing aid with the loop facility.

The user may then adjust the loop volume or tone to suit their own hearing preference without affecting others.



Personal Listeners

Personal Listeners help maintain conversation flow by amplifying the sounds you want to hear, and effectively reducing background noise.

With a personal listener, you can connect directly to various types of audio/ visual equipment, whether you want to use it during the day for conversation with friends or at work where loud noises might be too stressful. You can also use it to listen to your TV, radio, music device, or your computer.

You choose the way you want to listen, whether it's with a headset or neckloop, which you can use with your hearing aid via the "T" setting or loop program.

Crescendo 60/1

Under-chin stethoscope headset, plug-in microphone, microphone extension lead

Crescendo 60/1	
Order Code: 109-101	Colour: White/Grey
Price: (Ex. VAT) £79.99 Price: (Inc. VAT) £95.99	Size: 95 x 62 x 25mm
VAT Relief Yes	



Crescendo 60/5

Neckloop, plug-in microphone, microphone extension lead



Crescendo 60/5	
Order Code: 109-105	Colour: White/Grey
Price: (Ex. VAT) £79.99 Price: (Inc. VAT) £95.99	Size : 95 x 62 x 25mm
Price: (Inc. VAT) 293.99	

Alerting Signolux

' = ' / '

Relie

Signolux Alerting

Become aware of everyday occurances in the home or at work. E.g. ringing of the telephone, doorbell, fire alarm, or the sounds of a baby crying. "Feel" these occurances via a vibrating pad or "see" them happening via a flashing light.

Signolux Tower Receiver

Introducing this rechargeable Signolux Tower Receiver, which comes with rechargeable batteries and a power supply. Its aesthetic look means it can be positioned around the home & passed off as a decorative piece! Offering a modernised look, the buttons are larger and more obvious, which makes it easier to operate without having to pick it up, perfect for those who want a simple, fussfree experience.

Signolux Tower Receiver

Order Code: 230-713	Colour: White
Price: (Ex. VAT) £95.90 Price: (Inc. VAT) £115.08	Size: 143mm x 80mm

Signolux Doorbell

The Signolux Visual Alert system is a wireless, portable unit, perfect for use throughout the home (can be wall mounted). Operates up to 200 metres away from the Bell Push on the 868Mhz frequency. (*Portable/wall mounted receiver and bell push included*).

signelux

The receiver has 8 different indicators meaning that you can be alerted to different events within the home by adding different transmitters to the system. The receiver will alert you using very loud selectable chimes (>90dB), a bright flashing light and also visual indication symbols.

Volume control CD quality chime sound with choice of 15 sounds 8 visual notification icons White Strobe for visual indication



signelux

Signolux Phone Universal Transmitter

The Signolux Direct/Phone Universal Transmitter can be used to alert the user to the telephone ringing (by plugging into a BT socket), can be used to wire into an existing hard-wired doorbell system, alarm panel or water level indication.



Signolux	Phone	Transmitter

Order Code: 230-721	Colour: Black/Silver
Price: (Ex. VAT) £34.52 Price: (Inc. VAT) £41.42	Size : 122 x 75 x 38mm
VAT Relief Yes	

Signolux Audio Universal Transmitter

The Signolux Audio Universal Transmitter can be used to pick up audio signals from different sources such as a baby crying, an existing doorbell, an intercom system or the telephone ringing.

Signolux Audio	Universal Transmitter
Order Code: 230-720	Colour: Black/Silver
Price: (Ex. VAT) £34.52 Price: (Inc. VAT) £41.42	Size: 122 x 75 x 38mm
VAT Relief Yes	





Signolux Portable Receiver Pack

The Signolux Portable Receiver pack is made up of a Signolux portable vibrating receiver, Signolux charger cradle and Signolux vibrating pad.

Signolux Portable Receiver Pack

Order Code: 230-715	Colour: Black/Silver
Price: (EX. VAI) £135.90	Size: Portable Receiver: 90 x 58 x 30mm Charger: 75 x 80 x 65mm Vibrating Pad: 82 x 82 x 30mm



BABY MONITOR WITH VIBRATING PAD



AMPLICALL SENTINEL 1 VIDEO BABY MONITOR WITH LOUD ALARM AND VIBRATING PAD

NEW ARRIVAL! After much demand, we are proud to announce the arrival of this Amplicall Sentinel 1 Video Monitor, comprising of a baby camera unit, a parent screen receiver, and a wired vibrating pad.

The vibration feature and bright flashes act as an extra nudge for those with visual and hearing impairments. When the baby monitor is triggered, the receiver will ring, flash, vibrate and display the Live feed onto your highquality video receiver so you can sleep easy knowing your baby is too. Camera with night mode Motion detector Temperature detector Noise detector Feeding time reminder

Baby Monitor	
Order Code: 252-120	Colour: White
Price: (Ex. VAT) £148.99 Price: (Inc. VAT) £178.79	Size: 105 x 75 x 25mm
VAT Relief Yes	

Chimeflash Alerting



Chimeflash

alerting systems.

Range of up to 200m Easy to pair Wireless Installation

Chimeflash Doorbell

Extra loud and ultra-bright alerting, perfect for d/Deaf and hard of hearing people

This wireless, extremely loud doorbell alerting system offers ultra-bright flashing LED lights, volume control, & indicator symbols. Perfect for d/Deaf, hard of hearing, & visually impaired individuals, or anyone wanting a loud doorbell with excellent range & reliability!

Chimeflash Doorbell	
Order Code: 150-003	Colour: White
	Size: Transmitter: 80x 32 x 17mm Receiver: 145 x 105 x 52mm
VAT Relief Yes	

Chimeflash Sound Sensor

The new Chimeflash Sound Sensor will attach to any device which emits a loud repeating sound such as a doorbell, door entry system or alarm unit.

It will send a wireless signal to the Chimeflash Receiver (sold separately with our Chimeflash Doorbell), alerting you via extra loud sounds, ultra-bright flashing LEDs, and coloured indicators.

- Range of up to 200m
- Easy to pair
- Wireless Installation

Chimeflash Sound Sensor	
Order Code: 150-028	Colour: White
Price: (Ex. VAT) £27.65 Price: (Inc. VAT) £33.18	Size: 80 x 32 x 17mm
VAT Relief	

NEW

Chimeflash Telephone Transmitter

Range of up to 200m Easy to pair Wireless Installation

The Chimeflash Telephone Transmitter alerts you when your telephone is ringing via extra loud sounds, ultra-bright flashing LEDs and coloured indicators on your Chimeflash receiver.


You simply program this transmitter to the Chimeflash receiver via a couple of button presses and then plug the transmitter into the BT style phone socket using the double adapter provided.... now whenever the phone rings, so will your Chimeflash!

Chimeflash Telephone Transmitter

Order Code: 150-024	Colour: White
Price: (Ex. VAT) £27.65 Price: (Inc. VAT) £33.18	Size: 80 x 32 x 17mm

Alarm Clocks Alerting

Wake up on time every time with our extra loud alarm clocks, with flashing lights and a vibrating pad for that extra nudge.



Vintage Wake 'n' Shake

		IP
	9	1 a
	8	
	- 30	0 5
C		
0		

Vintage Wake 'n' Shake		
Order Code: 252-092	Colour: Black/White	
Price: (Ex. VAT) £49.99 Price: (Inc. VAT) £59.99	Size: 140 x 140 x 70mm	
VAT Relief Yes		

Wake 'n' Shake Curve

	Colour: Black
Price: (Ex. VAT) £44.99	Size: Clock: 178 x 94 x 50mm Pad: 90Ømm x 25mm
Price: (Ex. VAI) £44.99 Price: (Inc. VAT) £53.99	SIZE: Clock: 178 x 94 x 50 Pad: 90Ømm x 25mm



Wake 'n' Shake Curve



Phones Cordless

All of the phones in this section have either been specially designed for people with a hearing loss or have features that make them suitable. We have a selection of desktop (corded), and portable (cordless) phones, perfect for the home, office, or when out and about.



Xtra 3355 Combo

This telephone with answering machine is an easy-to-use amplified (+40dB) duo phone with a corded base and cordless DECT handset.

Xtra 3355 Combo	
Order Code: 255-627	Colour: White
Price: (Ex. VAT) £99.99 Price: (Inc. VAT) £119.99	Ringer : up to 90db
VAT Relief Yes	



2355 Xtra Handset

An extra handset for use or to be added to the Swissvoice range of Xtra cordless and combo phones in the 2355 & 3355 ranges.

2355 Xtra Handset		
Order Code: 255-629	Colour: White	
Price: (Ex. VAT) £33.33 Price: (Inc. VAT) £40.00	Talk Time - 10 Hours	
VAT Relief Yes		



Xtra 2355 Duo

The Xtra 2355 DUO is an easy-to-use amplified (up to 40dB) cordless twin handset phone with an answering machine and a call- block facility which is incredibly useful if you want to prevent nuisance calls.

Xtra 2355 Duo	
Order Code: 255-622	Colour: White
Price: (Ex. VAT) £91.66 Price: (Inc. VAT) £109.99	Talk Time - 10 Hours
Price. (Inc. VAT) £109.99	nours

AmpliDECT 295

This phone is hearing-aid compatible & has an adjustable volume of up to 30dB. Adjustable shaker function & a ringer volume provided. The AmpliDect 295 also comes supplied with a built-in answering machine that provides up to 15 minutes of recording time & remote access to your messages.

AmpliDECT 295	
Order Code: 252-040	Colour: White
Price: (Ex. VAT) £58.32 Price: (Inc. VAT) £69.98	Size : Handset: 172 x 55 x 25mm Base: 136 x 91 x 105mm
VAT Relief Yes 2 Vear	



AmpliDECT 595 Twin

This dual handset is ideal for anyone with hearing loss, dexterity, or visual impairments. This phone offers large backlit keys for easy dialling, receiving volume upto +50dB, ten one-touch number buttons for quick dial, & an SOS function, which offers added security to alert 4 of your contacts in an emergency.



AmpliDECT 595 Twin	
Order Code: 252-046	Colour: White
Price: (Ex. VAT) £133.32 Price: (Inc. VAT) £159.98	Size 72x55x25mm
VAT Relief Yes 2 ^{Vear}	

Phones Corded

AmpliPower 50

The AmpliPower 50 is the loudest and clearest amplified telephone available. With an amplification level of +60dB and the very useful hands-free speakerphone facility, tone, and memory functions, this telephone is suitable for anyone who finds telephone conversations difficult to hear.

AmpliPower 50	
Order Code: 252-195	Colour: White
Price: (Ex. VAT) £83.32 Price: (Inc. VAT) £99.98	Size : 195 x 210 x 70mm
VAT Relief Yes 2 ^{Year}	



Magna 4000

This phone is the elegant hands-free choice for anyone with moderate to severe hearing loss. It combines a very loud flashing ringer, large keys, caller ID, hands-free speakerphone, and more. Set the volume up to +60dB to get crystal clear sound quality.

Magna 4000	
Order Code: 253-026	Colour: Black
Price: (Ex. VAT) £70.82 Price: (Inc. VAT) £84.98	Ringer : up to 90db
VAT Relief Yes 2 ^{Year}	



CL100 Big Button Telephone

The CL100 is a corded telephone with a large button keypad and many useful features including the ability to amplify incoming speech (+30dB). The CL100 produces a loop signal which makes it suitable for use by people who wish to use their hearing-aid via the "T" setting.



Xtra 1110 Telephone

This easy-to-use amplified (+30dB) corded phone is made for people needing maximum comfort. Pin the 6 photo buttons with your relatives' pictures to instantly call them in the easiest way. In addition, the phone is fitted with 2 direct memory keys.

Xtra 1110 Telephone		
Order Code: 255-600	Colour: White	
Price: (Ex. VAT) £33.33 Price: (Inc. VAT) £40.00	Size : 217 x 200 x 82mm	
VAT Relief		



Tinnitus SnoozeBands

Compatible with many Bluetooth devices, these SnoozeBands connect to your mobile phone so you can listen to relaxing sounds to ease your tinnitus when tryng to sleep. The SnoozeBands range from 10-15-hour play time, offer soft padded speakers and a built-in microphone to answer calls.

SnoozeBand[™] Original Grey

Automatic shut-off timer after 1 or 2 hours, integrated microphone, and thin padded speakers

SnoozeBand ™	Original Grey
Order Code: 260-100	Colour: Grey
Price: (Inc. VAT) £38.99	Battery life: 8-10 hours from 1 charge
VAT Relief	









SnoozeBand[™] Pro 2.0

Blocks out 100% of light & includes built-in soft padded speakers

SnoozeBand™ Pro 2.0		
Order Code: 260-104	Colour: Black	
Price: (Inc. VAT) £44.99	Battery life: 8-10 hours from 1 charge	
VAT Relief No		

SnoozeBand[™] Deluxe

100% light-blocking & includes built-in microphone

SnoozeBand™	Deluxe	
Order Code: 260-107	Colour: Grey	
Price: (Inc. VAT) £49.99	Battery Life: 13-15 hours from 1 charge	
VAT Relief No		



How To Connect Your Devices to the TV

Below is a quick guide to help you set up your equipment with direct connection to Sky Q, Sky HD, & Virgin Media Set-Top Boxes via Optical Output. It covers the TV audio settings you should select and the audio outputs that are suitable for hearing devices like TV Listeners.

CONNECTING DIRECTLY TO A TELEVISION

If connecting directly to a TV set via optical then you must make sure that the TV is set to PCM. To do this go into the audio/sound settings menu on the TV and look for digital audio output/optical out. Normally this would be set to Dolby digital, it needs to be changed to PCM. Please note that some televisions auto switch this setting!

NOTE: Using the microphone supplied with your equipment attached to the speaker grille of the television will always return good results and is an excellent back-up way to connect to the TV.

SETTING UP YOUR EQUIPMENT WITH DIRECT CONNECTION TO SKY Q & SKY HD BOXES VIA OPTICAL

- Press Home on your Sky Q remote, select Settings, then select Setup
- Select Audio Visual, then select Digital Audio Output Optical and set it to 'NORMAL'
- Press the GREEN button to save the settings
- Power OFF the Sky Box, leave for 30 seconds and turn back on.
- Leave for a further 5 minutes until the box has had a chance to re-calibrate itself and then test the system

SETTING UP YOUR EQUIPMENT WITH DIRECT CONNECTION TO VIRGIN MEDIA SET TOP BOXES VIA OPTICAL

Press Home on your remote, select Help & Settings, then select Settings Select Audio > Dolby Digital to PCM. This setting applies to both HDMI and optical out. Power down the box for 5 minutes, restart the box, leave for a few minutes and test the system Virgin Media 360 Boxes ONLY: Select Settings > Audio & Video > PCM2.0 > Restart box as above.

WHY DON'T I GET AUDIO FROM DIGITAL SOURCES (NETFLIX, APPLE TV, FIRESTICK, ETC)?

Audio coming from certain content providers or devices (Netflix, Amazon Prime, Apple TV etc.) is set to **Output Dolby** and not **PCM** by default. To resolve this:

Change the 'Digital Audio Output' format in the app or device to PCM.

Then try the normal audio settings within the app (Netflix, Apple TV, FireStick, etc). Please bear in mind that sometimes these audio settings are only available once something is playing on screen.

PLEASE NOTE: This does **NOT** guarantee that you will hear the audio from these sources through your device, but it does give the **best** chance. **N.B.** When connecting directly to an external device such as SKY or Virgin Media box, you will not hear channels & smart services from the TV through the connected equipment.

Dolby Digital Audio Decoder Kit

A clever little device which will, in most cases, allow you to connect any of our Television Listening Systems (LA240, Sonumaxx & Earis) to any television, set top box, or audio output device without needing to change the settings on the television to output PCM signals. This is especially useful when the television or other audio source doesn't auto switch to PCM output or when you cannot find the appropriate settings in the menu systems.

Dolby Digital Audio Decoder Kit		
Colour: Black		
Size: 105 x 75 x 25mm		
(





Sarabec Arc Converter

If connecting to your television using our **ARC Audio extractor** then you must ensure that the HDMI socket you connect to is marked as **HDMI ARC or HDMI eARC**. If after connecting your equipment to the television via the ARC Audio extractor you do not get a sound output, or the sound you get is distorted then please change your television's audio output settings to **PCM, STEREO** or **AUTO** instead of Dolby Digital. Some televisions auto switch this function, most do not.

N.B Connection via the eARC socket will typically cut out the internal speakers from the television.



BRANDS WE WORK WITH:





Ordering Information

You can view and order all of our products online at www.sarabec.com 24 hours a day 7 days a week. We use a secure ordering system to ensure the safety of your details.

- If you are eligible for VAT relief, then, for each product required you only pay the price excluding VAT. You will need to complete the VAT relief certificate which can be found on our website.
- If a VAT exclusive price is not quoted for the product then the price including VAT must be paid.
- If paying by VISA, Mastercard, Maestro or AMEX please complete the credit/debit card payment details.
- Post or fax your order to us.

OR

• If paying by cheque, make it payable to 'Sarabec Ltd'. Enclose your payment with this order and post to us. Please note that cheques are cleared before goods are sent.

VAT

Throughout this catalogue you will find products that are priced inclusive (inc VAT) and exclusive (ex VAT) of VAT. VAT Relief is not available on postage and packing.

In the description of each product where we detail prices including and excluding VAT, that product is eligible for VAT relief. If an exclusive price is not shown, VAT relief is not available.

You are entitled to VAT relief if you are handicapped by a hearing disability. To benefit, simply complete the VAT Relief Certificate found on our website. If you are unable to make this declaration, you must include the VAT in your order.

Payment

You can pay by using a cheque or postal order payable to SARABEC Ltd or by using a credit card or debit card. If you wish to use a credit or debit card you may place your order and pay online via secure payment services or direct by telephone.

Trial Period

All products carry a 1-month trial period if purchased within the United Kingdom. You may return your products to us for any reason within 1 month of receipt and obtain a refund for your purchase provided that:

- You contact us first to let us know why you are returning the product.
- The products are returned in good condition and in their original packaging; we may charge if products have been misused.
- The products have been used in accordance with the product instructions.

The cost of returning any product is your responsibility. We cannot accept responsibility for loss or damage to goods until they arrive safely at our offices. We advise you to return goods adequately packed and covered by insurance where appropriate.

Delivery

Products are normally available from stock and are despatched within 48 hours of receipt of your order. Delivery to UK mainland addresses normally takes up to 3 working days from receipt of your order by us. Delivery is made through Royal Mail or independent carriers. A signature may be required on delivery of the product. Our delivery charges are for all UK mainland addresses, except the Channel Islands. For addresses outside this area please ask for delivery costs before ordering.

Further Information

Some products in this catalogue may not be available for the whole life of the catalogue.

All personal information entered on the Order or VAT relief forms will be treated as confidential as detailed in our Privacy Policy. You can ask for a copy of the Privacy Policy or obtain one from our website.

Our normal hours of business are: 8.15am to 4.45pm - Monday to Thursday, 8.15am to 4.00pm Friday.

Our full terms and conditions of sale are available on our website: www.sarabec.com.

All prices are liable to change (E&OE).

Sarabec Limited, 15 High Force Road, Middlesbrough, TS2 1RH, United Kingdom

Telephone (voice) 01642 247789 Next Generation Text (NGT) 18001 07497 242 294 Textphone (minicom) 01642 244148 Website www.sarabec.com e-mail mail@sarabec.co.uk



PLACE AN ORDER

To see the full range of our assistive listening equipment, or to place an order, visit: www.sarabec.com.



03/23

SARABEC LIMITED

15 High Force Road, Middlesbrough,TS2 1RH, UK TEL: 01642 247789 NGT: (Next Generation Text): 18001 07497 242 294 EMAIL: mail@sarabec.co.uk | WEBSITE: www.sarabec.com

JOIN US:

